

11-06-25

Yesterday each of us found a different Chl/Car, Chla/Chlb and chlorophyll concentration:

	Mads	Koen	Eva	Jasmya
Chl/Car	4.18	3.86	4.14	4.19
Chl a/Chlb	2.78	2.77	2.37	2.59
Chl concentration	0.63	1.59	0.86	1.24

↳ mg/mL

ABSORPTION

- Today we started with preparing our samples for an absorption spectrum. We prepared tubes with 3 mL B₁ and a bit (about half a milliliter) of our concentrated thylakoids.
- We moved to the spectrometer in the dark room and installed an integration sphere. (For measuring the scattering as well)
- We put our samples into the spectrometer and retrieved graphs we will analyse later today. We saved our samples for fluorescence.

FLUORESCENCE

- We determine fluorescence at 77 K (-196.15°C). The trapping to the RC will be less efficient at this temperature, so we will measure more fluorescence.
- To do this we retrieved liquid nitrogen from the -1 floor. We carried it up the stairs because turning it into the elevator is a safety hazard.
- We set up the machine for low temperatures.
- We froze our samples from the absorption in a thin glass pipette (we put them carefully into the liquid nitrogen, otherwise the glass would burst). Then we quickly transported the pipette into the machine (into a vacuum tube filled with liquid nitrogen).
- We retrieved our graphs and could directly see the fluorescence (between 730 nm and 950 nm). Data-analysis we will do this afternoon.

Data analysis

• we began making 8 graphs of our own plant and 3 graphs of all of our plants. In this way we can interpret our measurements. The graphs are normalized.

- ① Emission plant vs. absorption plant (room temp)
- ② Absorption plant vs absorption AT (model plant)
- ③ Emission plant vs emission AT (room temp) + difference
- ④ Emission plant vs emission AT (low temp) + difference
- ⑤ Emission plant room temp vs low temp
- ⑥ Absorption all plants
- ⑦ Emission all plants (room temp.)
- ⑧ Emission all plants (low temp.)

The code in python is uploaded in Github.

We didn't finish all the graphs today, so we finished them at home.

1.24 mg/mL, we need 0.5 mg

$$J \quad m = \frac{0.5}{1.24} = 0.4 \text{ mL}$$

$$K \quad \frac{0.5}{1.59} = 0.314 \text{ mL}$$

$$M \quad \frac{0.5}{0.63} = 0.79 \text{ mL}$$

$$E \quad \frac{0.5}{0.86} = 0.58 \text{ mL}$$

• We put these quantities into tiny tubes. pipet up and down to remix the thylacoids. Add 400 μ L of sd 1 ~~then~~

• 250 mL S3

• 52 to 125 mL (can be a bit more)

} See protocol
'Sucrose gradient'

We ~~finally~~ finished the day with putting the tubes in the ultra-centrifuge at 40000 RPM ($\sim 205000 g$) for 17h at 4°C. Tomorrow at 10:30 will continue with the gradient. One was a bit low concentrated and might fail. We discuss the graphs tomorrow morning.